


### NT-proBNP ELISA (Cat.No. SK-1204)

#### For the Determination of NT-proBNP (1-76) in Human Samples

#### **ASSAY CHARACTERISTICS**

<b>Method</b>	Sandwich ELISA, HRP/TMB, 12x8-well strips
<b>Antibodies/Standard</b>	Capture antibody: polyclonal anti human NT-proBNP antibody, specifically binding to amino acids 31-57 of proBNP. Detection antibody: polyclonal anti human NT-proBNP antibody, conjugated to peroxidase, specifically binding to amino acids 8-29 of proBNP. Standard material: synthetic human NT-proBNP (1-76).
<b>Sample type</b>	Human serum, EDTA plasma
<b>Standard range</b>	0-640 pmol/l (0, 10, 40, 160, 640)
<b>Conversion factor</b>	1 pmol/l = 8.475 pg/ml <i>refers to NT-proBNP (1-76) that is detected by the ELISA</i>
<b>Sample volume</b>	50 µl / well – direct measurement
<b>Sensitivity</b>	LOD (0pmol/l + 3SD): 3.0 pmol/l; LLOQ: 3.3 pmol/l
<b>Specificity</b>	This assay recognizes endogenous (natural) and recombinant human NT-proBNP (1-76).
<b>Incubation</b>	3 h /30 min – room temperature

#### **Typical standard curve of Biomedica NT-proBNP ELISA:**



**Serum values from apparently healthy donors:**

<b>Serum (n=70)</b>	<b>NT-proBNP (pmol/l)</b>
<b>Median</b>	<b>5.8</b>
Mean	9.8
Maximum	63.9
Minimum	1.1
Percentil 95%	32.5
Percentil 5%	2.0

**EDTA plasma values from apparently healthy donors:**

<b>EDTA plasma (n=28)</b>	<b>NT-proBNP (pmol/l)</b>
<b>Median</b>	<b>5.6</b>
Mean	9.5
Maximum	32.0
Minimum	0.6
Percentil 95%	27.9
Percentil 5%	0.8

**Serum values from unselected hospital panel:**

<b>Serum (n=117)</b>	<b>NT-proBNP (pmol/l)</b>
<b>Median</b>	<b>46.5</b>
Mean	84.5
Maximum	621.1
Minimum	0.0
Percentil 95%	22.1
Percentil 5%	0.8

**EDTA plasma values from unselected hospital panel:**

<b>EDTA plasma (n=40)</b>	<b>NT-proBNP (pmol/l)</b>
<b>Median</b>	<b>28.6</b>
Mean	57.9
Maximum	489.2
Minimum	0.3
Percentil 95%	157.7
Percentil 5%	1.2

**Serum values of a heart failure panel NYHA classification II-IV measured with the Biomedica NT-proBNP ELISA (#SK-1204)**

<b>Serum (n=63)</b>	<b>NT-proBNP (pmol/l)</b>
mean	54.0
max	399.3
min	1.1

**NYHA classification subgroups:**

measured with the Biomedica NT-proBNP ELISA (#SK-1204)

Data on serum values of a heart failure panel **NYHA classification II**

Serum (n=41)	NT-proBNP (pmol/l)
mean	25.5
max	75.6
min	1.1

Data on serum values of a heart failure panel **NYHA classification III**


Serum (n=14)	NT-proBNP (pmol/l)
mean	63.2
max	264.3
min	9.2

Data on serum values of a heart failure panel **NYHA classification IV**

Serum (n=7)	NT-proBNP (pmol/l)
mean	166.8
max	399.3
min	84.1

**MATRIX COMPARISON**

**Correlation of serum and plasma samples from apparently healthy individuals**



## **PERFORMANCE CHARACTERISTICS**

### **Spike Recovery:**

**The mean recovery of rec. NT-proBNP in serum samples is 99 and 108%.**

#### Experiment:

Recovery of spiked samples was tested by adding 2 concentrations of human recombinant NT-proBNP (1-76) to 4 different human serum samples.

Data showing spike/recovery of human serum samples:

<b>Matrix</b>		<b>Serum</b>			<b>S/R (%)</b>		
<b>Spike (pmol/l)</b>		<b>0</b>	<b>80</b>	<b>320</b>	<b>80</b>	<b>320</b>	
<b>Sample ID</b>	#1	1.6	88.5	347.2	109	108	
	#2	0.1	74.2	360.2	93	113	
	#3	1.7	80.9	338.7	99	105	
	#4	4.4	81.6	346.5	96	107	
					<b>Mean (%)</b>	<b>99</b>	<b>108</b>

**The mean recovery of rec. NT-proBNP in EDTA plasma samples is 93 and 94%.**

#### Experiment:

Recovery of spiked samples was tested by adding 2 concentrations of human recombinant NT-proBNP (1-76) to 4 different human EDTA plasma samples.

Data showing spike/recovery of human EDTA plasma samples:

<b>Matrix</b>		<b>EDTA plasma</b>			<b>S/R (%)</b>		
<b>Spike (pmol/l)</b>		<b>0</b>	<b>80</b>	<b>320</b>	<b>80</b>	<b>320</b>	
<b>Sample ID</b>	#1	0.8	73.4	276.8	92	87	
	#2	7.7	72.8	279.0	91	87	
	#3	9.3	77.2	324.1	97	101	
	#4	19.1	74.9	309.5	94	97	
					<b>Mean (%)</b>	<b>94</b>	<b>93</b>

### **Dilution Linearity:**

**Serum (n=3): 112-119%** (spiked with recombinant NT-proBNP)  
**79-86%** (endogenous NT-proBNP)

#### Experiment 1:

Dilution linearity was assessed by spiking serum samples with recombinant NT-proBNP (1-76) (320 pmol/l) and by diluting the samples with STD1 (serum matrix containing 0 pmol/l NT-proBNP, provided in the kit):

Data showing the dilution of recombinant NT-proBNP (1-76) in serum samples:

Serum Sample ID	Reference NT-proBNP (pmol/l)	Spiked NT-proBNP (pmol/l)	Dilution 1+1	
			NT-proBNP (pmol/l)	R (%)
#1	1.7	338.7	202.2	119
#2	0.1	360.2	202.4	112
#3	8.0	278.5	165.7	119
			<b>Mean (%)</b>	<b>117</b>

Experiment 2:

Dilution linearity was assessed by diluting serum samples containing endogenous NT-proBNP with STD1 (serum matrix containing 0 pmol/l NT-proBNP, provided in the kit).

Data showing the dilution of endogenous NT-proBNP (1-76) in serum samples:

Serum Sample ID	Reference NT-proBNP (pmol/l)	Dilution 1+1	
		NT-proBNP (pmol/l)	R (%)
#1	44.2	17.4	79
#2	51.3	20.4	80
#3	61.7	26.5	86
		<b>Mean (%)</b>	<b>82</b>

**EDTA plasma (n=2): 83-84%** (spiked with recombinant NT-proBNP)  
**80-85%** (endogenous NT-proBNP)

Experiment 1:

Dilution linearity was assessed by spiking EDTA plasma samples 1+1 with recombinant NT-proBNP (1-76) (640 pmol/l) and by diluting the samples with STD1 (serum matrix containing 0 pmol/l NT-proBNP, provided in the kit):

Data showing the dilution of recombinant NT-proBNP (1-76) in EDTA plasma samples:

EDTA plasma Sample ID	Reference NT-proBNP (pmol/l)	Spiked NT-proBNP (pmol/l)	Dilution 1+1	
			NT-proBNP (pmol/l)	R (%)
#1	36	347	145	84
#2	0.8	277	115	83
			<b>Mean (%)</b>	<b>84</b>

Experiment 2:

Dilution linearity was assessed by diluting EDTA plasma samples containing endogenous NT-proBNP with STD1 (serum matrix containing 0 pmol/l NT-proBNP, provided in the kit).

Data showing the dilution of endogenous NT-proBNP (1-76) in EDTA plasma samples:

	Reference	Dilution 1+1	
EDTA plasma Sample ID	NT-proBNP (pmol/l)	NT-proBNP (pmol/l)	R (%)
#1	269	114	85
#2	634	263	83
		<b>Mean (%)</b>	<b>84</b>

### **Intra-assay precision & Inter-assay precision:**

**Intra-assay (n=3) ≤ 4%, Inter-assay (n=8) ≤ 7%**

#### Experiment:

Intra-assay: 2 samples of known concentrations were tested 3 times in 1 assay by 1 operator.

Inter-assay: 2 samples of known concentrations were tested 8 times in 2 assays by different operators.

Data showing intra-assay and inter-assay precision:

Intra-assay (n=3)	Sample 1	Sample 2	Inter-assay (n=8)	Sample 1	Sample 2
Mean (pmol/l)	60.2	35.2	Mean (pmol/l)	52.1	108.1
SD (pmol/l)	2.0	0.9	SD (pmol/l)	1.7	7.9
CV (%)	4	3	CV (%)	3	7

### **The limit of quantification (LOQ):**

The LOQ is defined as the mean value of the back calculated concentration plus 3 times the standard deviation. The LOQ of the NT-proBNP ELISA is 3 pmol/l.

### **The lower limit of quantification (LLOQ):**

The lower limit of quantification is defined as the accuracy of the back calculated concentrations and shall not exceed ±25% (acc. to ICH (Ref. 1)).

For the NT-proBNP ELISA the LLOQ was determined as 3.3 pmol/l.

### **Specificity:**

This assay recognizes endogenous (natural) and recombinant human NT-proBNP (1-76).

### **Calibration:**

This immunoassay is calibrated against recombinant human NT-proBNP (1-76) peptide.

**SAMPLE CHARACTERISTICS**

**Sample stability:**

We recommend separating EDTA plasma or serum by centrifugation as soon as possible, e.g. 20 min at 2,000 x g, preferably at 4°C (2-8°C). EDTA-plasma or serum can be stored at 4°C (2-8°C) up to two days. For long term storage, aliquot the acquired EDTA plasma or serum samples and store them at -25°C or lower.

Samples can be subjected to 3 freeze-thaw cycles.


Serum and EDTA samples can be stored for at least 2 years at -80°C.

**Freeze/thaw of samples containing endogenous NT-proBNP**

The mean recovery of serum and EDTA plasma sample concentrations stressed by 3 F/T cycles is 94%. The mean CV of sample concentrations (not stressed and stressed up to 3 times by freeze-thaw cycles) is 4%.

NT-proBNP concentrations of samples after 3 freeze-thaw cycles:

Sample ID	NT-proBNP [pmol/l]				Mean	CV (%)	R (%) 3 F/T vs ref
	Ref	1	2	3			
Serum1	70	67	67	69	68	1	98
Serum2	15	13	15	14	14	5	98
Serum3	215	185	191	188	194	6	87
Serum4	141	129	129	126	131	4	90
Serum5	49	47	48	47	48	2	96
Plasma1	203	191	199	187	195	3	92
Plasma2	282	257	271	248	264	5	88
Plasma3	68	64	79	69	70	8	101
<b>Mean (%)</b>						<b>4</b>	<b>94</b>



### **Validation**

The assay is fully validated according to ICH Q2 (R1), Ref. (1).

### **References**

(1) CPMP/ICH/381/95 - ICH Topic Q2 (R1) „Validation of Analytical Procedures: Text and Methodology“ including:

ICH Q2A “Text on Validation of Analytical Procedures”

ICH Q2B “Validation of Analytical Procedures: Methodology”

### **Available on our homepage**

Enzyme immunoassay for the quantitative determination of NT-proBNP (1-76) in human serum.

Package Insert, MSDS, Information Folder, and References:

<http://www.bmgrp.com/products/cardiovascular/nt-probnp-elisa/>




### **NT-proBNP ASSAY COMPARISON**

#### **Comparison between NT-proBNP Assays: Biomedica SK-1204 versus Roche Elecsys 2<sup>nd</sup> generation Assay**


Data on serum values of a heart failure panel NYHA classification II-IV measured with the Biomedica NT-proBNP ELISA (#SK-1204):

n	63
mean	54.0 pmol/l
max	399.3 pmol/l
min	1.1 pmol/l




Data on serum values of a heart failure panel NYHA classification IV measured with the Biomedica NT-proBNP ELISA (#SK-1204)

n	7
mean	166.8 pmol/l
max	399.3 pmol/l
min	84.1 pmol/l



Data on serum values of a heart failure panel NYHA classification III measured with the Biomedica NT-proBNP ELISA (#SK-1204)

n	14
mean	63.2 pmol/l
max	264.3 pmol/l
min	9.2 pmol/l



**NT-proBNP ELISA, SK-1204 – ASSAY & PERFORMANCE CHARACTERISTICS**

Data on serum values of a heart failure panel NYHA classification II measured with the Biomedica NT-proBNP ELISA (#SK-1204)

n	41
mean	25,5 pmol/l
max	75,6 pmol/l
min	1,1 pmol/l

